## Restylane® Skinboosters™ Vital Light Lidocaine - Instructions for Use

**Composition**  
Hyaluronic acid, stabilized 12 mg/mL  
Lidocaine hydrochloride 3 mg/mL  
Phosphate buffered saline q.s.

**Description**  
Restylane Skinboosters Vital Light Lidocaine is a sterile, transparent gel of stabilized hyaluronic acid of non-animal origin with the addition of 0.3% lidocaine hydrochloride. It is supplied in a glass syringe. The product has a built-in dose-guide, Smart Click System, which when activated creates a clicking sound to indicate each injected dose. The 1 mL syringe gives approximately 100 doses. The contents of the syringe are sterilized using moist heat. The product is for single use only. Disposable 29G TW (thin-walled) needles, sterilized using ethylene oxide, are provided. To ensure traceability the patient record label (part of syringe label) should be attached to patient records.

**Intended use**  
This product is intended to improve the elasticity of the skin in the lower cheek/jawline in the face and upper neck. It should be injected in the dermal layer of the skin. The addition of lidocaine provides increased overall treatment comfort.

Before the first treatment session it is recommended to contact your local Galderma representative or Restylane distributor for more information about injection techniques and training opportunities. This product is only intended to be administered by authorized personnel in accordance with local legislation.

**Mode of action**  
This product is naturally integrated into the skin where it helps to improve the elasticity of the skin. This is accomplished by the water associated with the stabilized hyaluronic acid in the gel. The unique characteristics of the gel help maintain the effect for a long period of time.

**Performance**  
In clinical studies with Restylane Skinboosters Vital Light, patients experienced significant improvement in dermal elasticity after initial treatment.

**Contraindications**  
• Patients presenting with known hypersensitivity to hyaluronic acid filler, lidocaine or amide local anaesthetics.  
• Patients with a history of hypersensitivity to streptococcal proteins, as the product may contain trace amounts of such material.  
• Patients with bleeding disorders.  
• Patients presenting with porphyria.  
• Patients with known history of keloids.

**Warning**  
• Do not inject intravascularly. Introduction of product into the vasculature may lead to embolization, occlusion of the vessels, ischemia, or infarction.  
• Rare but serious adverse events associated with the intravascular injection of soft tissue fillers in the face have been reported and include temporary or permanent vision impairment, blindness, cerebral ischemia or cerebral hemorrhage, leading to stroke, skin necrosis, and damage to underlying facial structures.

• Immediately stop the injection if a patient exhibits any of the following symptoms, including changes in vision, signs of a stroke, blanching of the skin, or unusual pain during or shortly after the procedure.  
• Patients should receive prompt medical attention and possibly evaluation by an appropriate health care practitioner specialist should an intravascular injection occur.

• Defer use of the device where there is active disease such as inflammation, infection or tumours, in or near the intended treatment site until the disease has been controlled.  
• Do not resterilize.  
• Do not mix with other products prior to injection of the device.

### Precautions

**General considerations relevant to injectable medical devices**  
• In order to minimize the risks of potential complications, this product should only be used by health care practitioners who have appropriate training, experience, and who are knowledgeable about the anatomy at and around the site of injection. Knowledge of the anatomy of treatment site and special caution are required in order to avoid perforation or compression of vessels, nerves and other vulnerable structures.

• Health care practitioners are encouraged to discuss all potential risks of soft tissue injection with their patients prior to treatment and ensure that patients are aware of signs and symptoms of potential complications.  
• Injection procedures are associated with a risk of infection. Aseptic technique and standard practice to prevent cross-infections are to be observed.

• Use with caution in patients who are immunosuppressed. Special caution should be exercised when treating areas in close proximity to permanent implant.  
• Localized ischemia/necrosis and scarring may occur after injection in or near vessels. Special caution should be taken if the patient has undergone a prior surgical procedure in the planned treatment area.

• Special care should be exercised to avoid a too large or too superficial injection in treating facial areas with limited soft tissue support or soft tissue cover, to avoid formation of palpable lumps.  
• Injection procedures can lead to reactivation of latent or subclinical herpes viral infections.

• Bruising and bleeding may occur at injection sites. The device should be used with caution in patients who have undergone therapy with thrombolytics, anticoagulants or inhibitors of platelet aggregation, such as aspirin or non-steroidal anti-inflammatory drugs, in the preceding 2 weeks.

• Individuals with skin types (IV-VI) may be prone to keloid formation. Ensure that the patient has no history of previous keloid formation.  
• Patients with unattainable expectations are not suitable for treatment with this product.

• Do not use the product if package is damaged.  
**Specific considerations relevant to the use of this product**  
• Do not inject this product into an area where an implant of non hyaluronic acid has been placed.  
• Patients should avoid excessive sun or extreme cold at least until any initial swelling and redness has resolved.

• If laser treatment, chemical peeling or any other procedure based on active dermal response is performed after treatment with this product there is a theoretical risk of eliciting an inflammatory reaction at the injection site. This also applies if the product is administered before the skin has healed completely after such a procedure.

• Temporary palpable lumpiness has been noted in some patients. Do not inject too large volumes or too superficially, as the product then may cause intradermal lumps.  
• This product has not been tested in pregnant or breastfeeding women or in children.

• Considerations should be given to the total dose of lidocaine administered if dental block or topical administration of lidocaine is used concurrently. High doses of lidocaine (more than 400 mg) can cause acute toxic reactions manifesting as symptoms affecting the central nervous system and cardiac conduction.

• Lidocaine should be used with caution in patients receiving other local anaesthetics or agents structurally related to amide-type local anaesthetics e.g., certain anti-arrhythmics, since the systemic toxic effects can be additive.

• Lidocaine should be used cautiously in patients with epilepsy, impaired cardiac conduction, severely impaired hepatic function or severe renal dysfunction.

• Peribulbar injections of local anaesthetics carry a low risk of persistent ocular muscle dysfunction.

### Adverse events

**Anticipated injection-related reactions**  
Injection-related reactions might occur. These reactions include: bruising, erythema, itching, swelling, pain or tenderness at the implant site. Typically resolution is spontaneous within a few days after injection into the skin.

**Post marketing adverse event reporting**  
The frequency of post market adverse event reporting is calculated on the estimated number of treatments performed with the Restylane Skinboosters range of products.

I/10 000 – I/50 000: Erythema, inflammation, pain/tenderness, papules/nodules, swelling  
I/50 000 – I/100 000: Bruising, induration  
<I/100 000: Abscess, acne, atrophy/kyriarring, urticaria, dermatitis, discoloration, granuloma, blisteritis, infection, ischemia/necrosis, mass, neurological symptoms such as paresthesia, pruritus, rash, reactivation of herpes infection, short duration of effects, telangiectasia, urticaria

• Vascular compromise may occur due to an unintentional intravascular injection or as a result of vascular compression associated with implantation of any injectable product. This may manifest as blanching/discoloration such as a dusky or reticular appearance of the tissue, necrosis or ulceration at the implant site or in the area supplied by the blood vessels affected; or rarely as ischemic events in other organs due to embolization. Rare but serious cases of ischemic events associated with temporary or permanent vision impairment, blindness, cerebral ischemia or stroke have been reported following facial aesthetic treatments.

In rare cases intradermal lumps have been reported to remain for several months or very rarely longer than one year.  
As for injectable dermal fillers, symptoms of inflammation at the implant site commencing either shortly after injection or after a delay of 2-4 weeks have been reported. In case of unexplained inflammatory reactions infections should be excluded and treated if necessary (antibiotics inadequately treated infections may progress into complications such as abscess formation). Treatment using only oral corticosteroids without concurrent antibiotic treatment is not recommended.

The prolonged use of any medication, e.g. corticosteroids or antibiotics in treatment of adverse events has to be carefully assessed, since this may carry a risk for the patient. In case of persistent or recurrent inflammatory symptoms, consider removal of the product by aspiration/drainage, extrusion or enzymatic degradation (hyaluronidase injections). Before any removal procedure is performed, the swelling may be reduced by using e.g. NSAID for 2-7 days or a short course of corticosteroids for less than 7 days, in order to more easily palpate any remaining product.

For patients who have experienced clinically significant reactions, a decision for retreatment should take into consideration the cause and significance of previous reactions.

Post inflammatory pigmentation changes have been observed in clinical studies in people with dark skin (Fitzpatrick Type IV-VI).  
For reporting of adverse event contact your local Galderma representative or Restylane distributor for this product.

### Needle

Disposable sterile 29G TW needles are provided. Alternatively, a sterile blunt cannula 30G can be used. The size and the length of the cannula will affect the force needed to extrude the gel. If a thinner cannula is used the resistance during injection may be too high resulting in an increased risk for leakage or separation of the cannula and syringe. The same considerations are applicable for needles.

### Assembly of the syringe

Use the thumb and forefinger to hold firmly around both the glass syringe and the luer-lock adapter. Grasp the needle shield (or hub if using cannula) with the other hand. To facilitate proper fastening, both **push and rotate** firmly, see Figure 1. Strict aseptic technique must be followed. Improper assembly may result in separation of the needle and syringe during injection.

### Treatment procedure

• The patient shall be informed about the indications, expected result, precautions and potential adverse events. The patient's need for additional pain relief should be assessed.  
• Clean the treatment site thoroughly with a suitable antiseptic solution.

• To avoid breakage of the needle or cannula do not attempt to bend or otherwise manipulate it before or during treatment.  
• This product should be injected into mid dermis of the skin.

• The Smart Click System is activated by pressing down the button located on the finger grip until it locks into place, see Figure 2.

• Before injecting, remove the air by pressing the rod carefully until a small droplet is visible at the tip of the needle. When the Smart Click System is switched on, press the plunger rod carefully until the first click is heard to prime the system before use.

• As an alternative to the needle, a blunt cannula can be used. After preparation as described above, an entry point is made in the skin, e.g. with a sharp needle of appropriate size. Inject slowly. During injection, keep the side hole of the cannula facing downwards, away from the skin surface, to ensure that the flow of the gel is maintained at the correct tissue depth.

• Inject the product slowly and apply the least amount of pressure necessary.  
• Do not apply excessive pressure to the syringe at any time. Presence of scar tissue may impede advancement of the cannula/needle. If resistance is encountered the cannula/needle should be partially withdrawn and repositioned or fully withdrawn and checked for function.

• If the overlying skin turns a whitish colour (blanching), the injection should be stopped immediately and the area massaged until it returns to a normal colour. Blanching may represent a vessel occlusion. If normal skin coloring does not return, do not continue with the injection.  
• A too large volume or a too superficial injection may give bumps on the treatment site.

• Treated areas can be gently massaged immediately after the injection if any irregularities are noted.  
• A treatment plan for this product with three treatments 2-4 weeks apart is recommended. Generally, a maintenance treatment is repeated every 4-6 months, but results and patient preferences may vary. The maximum volume of gel injected should not exceed 3.5 mL per treatment session and 17.5 mL per year and patient.

The syringe with finger grip and plunger rod, disposable needle/blunt cannula and any unused material must be discarded immediately after the treatment session and must not be reused due to risk for contamination of the unused material and the associated risks including infections. Disposal should be in accordance with accepted medical practice and applicable national, local or institutional guidelines.

### Shelf life and storage

The expiry date is indicated on package. Store up to 25 °C. Protect from freezing and sunlight.

### Manufacturer

Q-Med AB, Seminariegatan 21, SE-752 28 Uppsala, Sweden  
Phone: +46(0)18 474 90 00. Fax: +46(0)18 474 90 01  
www.q-med.com, e-mail: info.q-med@galderma.com


### Manufactured for

Galderma Canada Inc.  
Thornhill, ON L3T 7V9  
1 800 467-2081

Restylane and Galderma are trademarks of Nestlé Skin Health S.A.

**Figure 1.**  
Use the thumb and forefinger to hold firmly around both the glass syringe barrel and the luer-lock adapter. Grasp the needle shield (or hub if using cannula) with the other hand. To facilitate proper assembly, both **push and rotate** firmly.

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**Figure 2.**  
Activate the Smart Click System by pressing down the button on the finger grip (Figure 2a) until it locks into place (Figure 2b).

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Q-Med AB, Seminariegatan 21, SE-752 28 Uppsala, Suède
Restylane et Galderma sont des marques de commerce de Nestlé Skin Health S.A.