**RESTYLANE FINE LINES-L™ INSTRUCTIONS FOR USE**

**Formula**  
 Contient de l'acide hyaluronique stabilisé 20 mg/ml  
 Hyaluronic acid, stabilized 20 mg/ml  
 Lidocaïne hydrochloride 0,3 %  
 Phys. sodium chloride solution, pH 7 q.s.

**Description**  
*Restylane Fine Lines-L* is a sterile, transparent biodegradable gel supplied in a glass syringe together with a 30G or 29G Thin Wall (TW) needle. The product is for single use only. *Restylane Fine Lines-L* is a unique form of non-animal, stabilized hyaluronic acid (NASHA™). Hyaluronic acid is a natural polysaccharide which occurs as an important structural element in the skin and in subcutaneous and connective tissues as well as in the synovial tissue and fluid. Hyaluronic acid belongs to a group of very few substances which are identical in all living organisms.

**Modes of action**  
*Restylane Fine Lines-L* acts by adding volume to the tissue, thereby restoring the skin contours to the desired level of correction. *Restylane Fine Lines-L* is naturally integrated into the tissue and will in time undergo isovolumic degradation.

*Restylane Fine Lines-L* has been enhanced with the addition of lidocaine in order to reduce patient discomfort during treatment.

**Indication and usage**  
*Restylane Fine Lines-L* is intended to be used for facial tissue augmentation. It is recommended to use the product for the correction of thin superficial lines, for example periorbital lines, perioral lines, forehead lines and smile lines. It should be injected in the superficial dermal layer of the facial skin. Deeper injections into the subcutaneous fatty tissue or supraperiosteal layer are appropriate for areas with adequate soft tissue support and soft tissue cover such as midface and jaw line. With cutaneous contour deformities the best results are obtained if the defect can be manually stretched to the point where it is eliminated. The degree and duration of the correction depend on the character of the defect treated, the tissue stress at the implant site, the depth of the implant in the tissue and the injection technique. Markedly indurated defects may be difficult to correct. For the correction of cutaneous contour deformities such as wrinkles it is recommended that *Restylane Lidocaïne* be used. *Restylane Lyft Lidocaïne* can be used for shaping the contours of the face and for the correction of folds. For lip enhancement both *Restylane Lidocaïne* and *Restylane Lyft Lidocaïne* can be used. Also combinations of *Restylane Lidocaïne*, *Restylane Lyft Lidocaïne* and *Restylane Fine Lines-L* can be used. Please consult the *Restylane Lidocaïne* and *Restylane Lyft Lidocaïne* Instructions for Use for more information.

**Contraindications**  
 Patients with known hypersensitivity to hyaluronic acid filler, lidocaine or amide type local anaesthetics. Patients with bleeding disorders.

- Warnings**
- Do not inject intramuscularly or intravascularly. Introduction of product into the vasculature may lead to embolization, occlusion of the vessels, ischemia, or infarction. Aspiration prior to injection is recommended.
  - Rare but serious adverse events associated with the intravascular injection of soft tissue fillers in the face have been reported and include temporary or permanent vision impairment, blindness, cerebral ischemia or cerebral hemorrhage, leading to stroke, skin necrosis, and damage to underlying facial structures.
  - Immediately stop the injection if a patient exhibits any of the following symptoms, including changes in vision, signs of a stroke, blanching of the skin, or unusual pain during or shortly after the procedure.
  - Patients should receive prompt medical attention and possibly evaluation by an appropriate health care practitioner specialist should an intravascular injection occur.
  - Do not use in patients who are taking anticoagulants or antiplatelet drugs such as aspirin, NSAID for 2-7 days or a short course of corticosteroids for less than 7 days, in order to more easily palpate any remaining product.
  - Do not use in patients with a history of hypersensitivity to streptococcal proteins, as the product may contain trace amounts of such material.
  - Do not sterilize.
  - Do not mix with other products prior to injection of the device.
  - Do not use where there is active disease, such as inflammation, infection or tumours, in or near the intended treatment site.

**Precautions**  
 General considerations relevant to injectable medical devices

- In order to minimize the risks of potential complications, this product should only be used by health care practitioners who have appropriate training, experience, and who are knowledgeable about the anatomy at and around the site of injection.
- Health care practitioners are encouraged to discuss all potential risks of soft tissue injection with their patients prior to treatment and ensure that patients are aware of signs and symptoms of potential complications.
- Injection procedures are associated with a risk of infection. Aseptic technique and standard practice to prevent cross-infections are to be observed.
- Use with caution in patients who are immunosuppressed.
- Special caution should be exercised when treating areas in close proximity to permanent implants.
- Special caution is required to avoid perforation or compression of vessels, nerves and other vulnerable structure.
- Localized ischemia/necrosis and scarring may occur after injection in or near vessels, especially in areas with limited collateral circulation. Special caution should be taken if the patient has undergone a prior surgical procedure in the planned treatment area.
- Special caution should be exercised in treating facial areas with limited soft tissue support or soft tissue cover, such as the periorbital area, to avoid formation of palpable lumps.
- Patients with pre-existing pigmented dark lower eye lid circles, thin skin and pre-existing tendency toward edema formation are not suitable candidates for treatment of the lower periorbital region.
- Injection procedures can lead to reactivation of latent or subclinical herpes viral infections.
- Patients who are using substances that affect platelet function, such as aspirin and non-steroidal anti-inflammatory drugs may, as with any injection, experience increased bruising or bleeding at injection sites.
- Patients with unattainable expectations are not suitable candidates for treatment.
- Do not use the product if package is damaged. Specific considerations relevant to the use of this product
- Do not inject this product into an area where an implant other than hyaluronic acid has been placed.
- If the product is injected too superficially this may result in visible lumps and/or bluish discoloration. This may include avoid excessive sun or extreme cold at least until any initial swelling and redness has resolved.
- If laser treatment, chemical peeling or any other procedure based on active dermal response is performed after treatment with this product there is a theoretical risk of eliciting an inflammatory reaction at the implant site. This also applies if the product is administered before the skin has healed completely after such a procedure.
- This product has not been tested in pregnant or breastfeeding women or in children.
- Considerations should be given to the total dose of lidocaine administered if dental block or topical administration of lidocaine is used concurrently. High doses of lidocaine (more than 400mg) can cause acute toxic reactions manifesting as symptoms affecting the central nervous system and cardiac conduction.
- Lidocaine should be used with caution in patients receiving other local anaesthetics or agents structurally related to amide-type local anaesthetics e.g., certain anti-arrhythmic, since the systemic toxic effects can be additive.
- Lidocaine should be used cautiously in patients with epilepsy, impaired cardiac conduction, severely impaired hepatic function or severe renal dysfunction.
- Periubaric injections of local anaesthetics carry a low risk of persistent ocular muscle dysfunction.

**Adverse Events**  
 Anticipated injection-related reactions  
 Injection-related reactions might occur. These

reactions include bruising, erythema, itching, swelling, pain, or tenderness at the injection site. Typically resolution is spontaneous within a few days after injection into the skin.

**Post marketing adverse event reporting**  
 The following post marketing adverse events have been reported (non-exhaustive list). The frequency of reporting is based on the number of estimated treatments performed with *Restylane Fine Lines* products.

1/1 000 - 1/10 000	Swelling with immediate onset and onset up to several weeks after treatment
1/10 000 - 1/100 000	Mass/induration, Erythema, Pain/tenderness, Bruising/hematoma, Papules/nodules, Inflammation, Ischemia/necrosis, Short duration of effect, Infection and abscess formation, Discoloration, Hypersensitivity, Neurological symptoms including paraesthesia.
<1/100 000	Device dislocation, Injection site reactions including burning sensation, Pruritus, Atrophy/scarring, Capillary disorders including telangiectasia, Eye pain, Granuloma

Other potential adverse events that have been reported following injection of other *Restylane* products and that may occur when using *Restylane Fine Lines-L* include the following: Acne, Blisters/vesicles, Dermatitis, Fistula/dischage, Rash, Reactivation of herpes infection, Urticaria, and Visual disturbance including blindness, transient blurred vision, and reduced visual acuity. Post inflammatory pigmentation changes have been observed in clinical studies on other *Restylane* products in people with dark skin (Fitzpatrick Type IV - VI).

Vascular compromise may occur due to an unintentional intravascular injection or as a result of vascular compression associated with implantation of any soft tissue filler in the face. This may manifest as blanching, discoloration such as dusky or reticular appearance of the tissue, necrosis or ulceration at the implant site or in the area supplied by the blood vessels affected; or rarely as ischemic events in other organs due to embolization. Rare but serious cases of ischemic events associated with temporary or permanent vision impairment, blindness, cerebral ischemia or stroke have been reported following facial aesthetic treatments.

Localized ischemia/necrosis and scarring may occur after injection in or near vessels, such as in the lips, glabellar area, or in the nose. For patients who had prior rhinoplasty special caution should be taken in the nose area.

Symptoms of inflammation at the implant site commencing either shortly after injection or after a delay of up to several weeks have been reported. In case of unexplained inflammatory reactions infections should be excluded and treated if necessary since inadequately treated infections may progress into complications such as abscess formation. Treatment using only oral corticosteroids without concurrent antibiotic treatment is not recommended.


The prolonged use of any medication, e.g., corticosteroids or antibiotics in treatment of adverse events has to be carefully assessed, since this may carry a risk for the patient. In case of persistent or recurrent inflammatory symptoms, consider removal of the product by aspiration/drainage. Extrusion or enzymatic degradation use of hyaluronidase has been described in scientific publications. Before any removal procedure is performed, the swelling may be reduced by using, e.g., NSAID for 2-7 days or a short course of corticosteroids for less than 7 days, in order to more easily palpate any remaining product.

For patients who have experienced clinically significant reactions a decision for retreatment should take into consideration the cause and significance of previous reactions.


Adverse events must be reported to Galderma Canada Inc.

**Interactions**  
 Treatment with *Restylane Fine Lines-L* in combination with drugs and other devices or concomitant dermal therapies has not been evaluated in controlled clinical studies.

**Assembly of 30G needle to syringe**  
 For safe use of *Restylane Fine Lines-L* it is important that the needle is properly assembled. See pictures A and B.

**Assembly of 29G (TW) needle to syringe**  
 Use the thumb and forefinger to hold firmly around both the glass syringe barrel and the Luer-Lok adapter. Grasp the needle shield with the other hand. To facilitate proper assembly, both push and rotate firmly.

**Dosage and administration**  
 Before the treatment, the patient's suitability for the treatment and the need for pain relief should be assessed. Normally, no anaesthesia is necessary when treating thin superficial lines. The patient should be informed about the indications, expected result, contraindications, precautions, warnings and potential adverse events. The treatment site should be cleaned with an astringent antiseptic solution. *Restylane Fine Lines-L* is administered using a 29G (TW) or 30G needle by injecting the material into the superficial layer of the skin. An injection too superficial may give blanching effects and bumps on the treatment site. Aspiration prior to injection is recommended to verify that the needle is not intravascular. If the overlaying skin turns a whitish colour (blanching), the injection should be stopped immediately and the area massaged until it returns to a normal colour. Blanching may represent a vessel occlusion. If normal skin coloring does not return, do not continue with the

injection. Before injecting, press the rod carefully until a small droplet is visible at the tip of the needle. In order to minimize the risk of potential complications, inject the product slowly and apply the least amount of pressure necessary.

The injection technique with regard to the depth of injection and the administered quantity may vary. The linear threading technique can be used to carefully lift up the thin superficial line, but some physicians prefer a series of punctual injections or a combination of the two. During injection it is recommended that the eye of the needle should face upwards. Inject *Restylane Fine Lines-L* while pulling the needle slowly backwards. Injection should stop just before the needle is pulled out from the skin to prevent material from leaking out from the injection site. Please consult Galderma Canada Inc. for details. Defects should be fully corrected, but not overcorrected, at each treatment session. If the skin of the patient is very loose, it is recommended that *Restylane Fine Lines-L* is injected on two separate occasions. The correction site should be massaged to conform to the contour of the surrounding tissues. For each treatment site a maximum dosage of 0.8 ml is recommended. If the treated area is swollen directly after the injection, melting ice or a cold pack can be applied on the site for a short period. After the first treatment, additional implantations of *Restylane Fine Lines-L* may be necessary to achieve the desired level of correction. Periodic follow-up injections help sustain the desired degree of correction.

**Note!** The correct injection technique is crucial for the final result of the treatment. *Restylane Fine Lines-L* is only intended to be administered by authorized personnel in accordance with local legislation.

**The syringe, the needle and any unused material must be discarded directly after the treatment session.**

**How supplied**  
*Restylane Fine Lines-L* is supplied in a disposable glass syringe. *Restylane Fine Lines-L* is co-packed with a sterilized needle(s) as indicated on the carton, either 30 G x 1/2" or 29 G (TW) x 1/2". A patient record label is a part of the syringe label (see picture C). Remove it by pulling the flap marked with three small arrows (see picture D). This label is to be attached to patient records to ensure traceability of the product. The contents of the syringe are sterile. The number of units per package and the volume contained in each syringe is as stated on the outer package.

**STERILE NEEDLES**


- Follow national, local or institutional guidelines for use and disposal of medical sharp devices. Obtain prompt medical attention if injury occurs.
- Do not avoid breakage of the needle, do not attempt to bend or otherwise manipulate it before or during treatment. Discard it and complete the procedure with a replacement needle.
- Do not re-sterilise used needles. Recapping by hand is a hazardous practice and should be avoided.
- Discard unshielded needles in approved sharps collectors.

**Shelf life and storage**  
 As indicated on package. Store up to 25°C (77°F). Protect from sunlight and freezing. Refrigeration not required. Canadian Patent No. 2.226.488

**Manufactured by**  
 Q-Med AB  
 Seminariegatan 21, SE-752 28 Uppsala, Sweden  
 Phone: +46(0)18 474 90 00  
 Fax: +46(0)18 474 90 01  
 www.q-med.com  
 Email: info.q-med@galderma.com

**Manufactured for**  
 Galderma Canada Ltd.  
 Thornhill (Ontario) L3T 7V9  
 1 800 467-2081

**IF THE PACKAGE IS DAMAGED, DO NOT USE**

*Restylane* and *Galderma* are trademarks of Nestlé Skin Health S.A.

All other trademarks are the property of their respective owners.

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